



VERTICAL MINI DUAL GEL SYSTEM

BIO 05-02

Complete Solution for Life Science Research

VERTICAL MINI DUAL GEL SYSTEM (BIO 05-02)



Technical Specifications

Unit Dimensions (L x W x H)	16 x 20 x 15cm
Gel Dimensions (L x W)	8 x 7 cm (2 Gel)
Buffer Volume	150 + 200 ml
Maximum Sample Capacity	14
Combs	7 Well Teflon Comb 0.5mm,1mm – 2 No's Each
Spacers	Teflon Spacers 0.5mm,1mm – 4 No's Each
Platinum Electrodes	Red and Black (Fixed)
Glass Plates	Notched & Rectangular (2 Sets)
Gasket	Fixed
Clamps and Screws	2 Sets
Connecting Cord	1 Pair
Water Circulation	Yes
Unit Dimensions (L x W x H)	16 x 20 x 15cm
User Manual	

BIO WORLD

1247, 6th Main, 5th Stage, 1st Phase, B.E.M.L Layout , R.R Nagar, Bangalore – 560098

Website: www.bioworldproducts.com ; Email : info@bioworldindia.in

Mobile: +91-9591779325,7760658542,9448351365