



DRY BATH

BIO DB-01

Complete Solution for Life Science Research

DRY BATH



Dry Bath is Specially designed to cater the needs of laboratory requirements. The product is equipped with some user friendly like compact style, temperature control and convenient operational modes. The Company manufactures them in fulfill the needs of the customers who are associated with the laboratory testing application. We provide dry bath at a competitive market price. They are excellent to uses in the research laboratories. Dry Bath is widely used Biotechnology Instrumentation.

Specialty:

- High Quality
- Great Performance
- Compact Style
- Reliable



Technical Specification

Description	Dry Bath
Cabinet	MS Powder Coated
Block Size	9.5 X 7.5 X 5 Cms for 24 places
Anodized Aluminum block	8 X 1.5 ml, 8 X 0.5 ml & 8 X 0.2 ml Eppendorf tubes
Temperature	Ambient to 99° C/ $\pm 1^{\circ}$ C accuracy (PID Control)
Block	Anodized Aluminum Block
Sensor	PT 100 Sensor
Control	Electronic Controls with PID Display
User Manual	Yes

Ordering Information:

Cat#	Description	Qty
BIO DB-01	Dry Bath with Standard Aluminum Block	1 No

BIO WORLD

1247, 6th Main, 5th Stage, 1st Phase, B.E.M.L Layout , R.R Nagar, Bangalore – 560098

Website: www.bioworldproducts.com ; Email: info@bioworldindia.in

Mobile: +91-9591779325,7760658542,9448351365