



INCUBATOR ORBITAL SHAKER COOLING

Complete Solution for Life Science Research

Incubator Orbital Shaker Cooling



Incubator shaker is equipment designed with heating and cooling facility inside the chamber with orbital shaking platform. Some applications require shaking facility in a closed chamber with a set value of uniform temperature. Soft Switch controller with LCD display for temperature & RPM Display

Orbital incubator shaker cooling comes with different platform sizes, customer can choose as per there requirement.

Technical Specification

Description	Incubator Orbital Shaker Cooling
Outer Chamber	MS Powder Coated
Inner Chamber	Stainless Steel
Platform	Universal Spring Clamp
Controller	Soft Switch
Display	LCD Display
Movement	Orbital Motion
Temperature	10°C to 60°C
Speed	40 to 250 Rpm
Operating Voltage	230v ac@ 50 Hrtz
User Manual	Yes

Ordering Information

Cat#	Description	Qty
BIO-IOS-04 (a)	Incubator Orbital Shaker Cooling 25 x 25 cm	1 No
BIO-IOS-04 (b)	Incubator Orbital Shaker Cooling 38 x 28 cm	1 No
BIO-IOS-04 (c)	Incubator Orbital Shaker Cooling 45 x 45 cm	1 No

BIO WORLD

1247, 6th Main, 5th Stage, 1st Phase, B.E.M.L Layout , R.R Nagar, Bangalore – 560098

Website: www.bioworldproducts.com ; Email : info@bioworldindia.in

Mobile: +91-9591779325,7760658542,9448351365